

# Application

**Determination of isocyanate  
(NCO-) content**

# Application

## Use

This method is applicable for material containing reactive isocyanate groups. The isocyanate is reacted with di-n-butyl amine to form a urea. Unreacted (excess) amine is determined by backtitration with hydrochloric acid.

## Appliances

- Titrator: TL 7000/TL 7750 M1/20 or 50
- Basic device
- Magnetic stirrer TM 235
- 20/50 mL exchange unit WA 20/50, with amber glass bottle for the titrant, complete
- Erlenmeyer flask 250 or 300 ml, Syringe 2 or 5 ml for high viscosity samples

## Electrodes

- Electrode: N 6480 eth + TZ 1643 titration tip
- Electrolyte: L 5034 (LiCl/ethanol)
- Calibration: n.a.

## Reagents

- Titrant: HCl 1 mol/l aqueous or HCl in methanol. Sometimes HCl 0.5 mol/l
- Titer determination: TRIS (Tris hydroxyl aminomethan)
- Reagent mixture: Di-n-butylamine solution 1 mol/l
- Solvent: Toluene dried with molecular sieve 4A
- Additional solvent: acetone or other suitable solvents

## Description

### Determination of the exact concentration of the HCL titrant

We recommend ready to use HCl titrant. The exact concentration of the HCL 1.000 mol/l can be determined using the titrimetric standard Tris hydroxy amino methane (TRIS). The standard can be dried for 24 h in an exsiccator at room temperature.

1-2 g of the standard are weighed accurately in a 150 mL beaker and dissolved in 80 mL of dist. water with stirring. It is titrated with the 1 mol/l HCl solution.

Repeat the standardization two times. The average value is stored automatically in the exchangeable unit.


## Application

Page 1: Curve and result: Titer determination

The curve shows the use of an HCl 0.1 mol/l. The sample amount of TRIS should be 10 – 15 times higher for the 1 mol/l titrant.

### GLP documentation

#### Titration graph



#### Method data

|              |           |                     |          |
|--------------|-----------|---------------------|----------|
| Method name: | Titre HCl | Titration duration: | 3 m 8 s  |
| End date:    | 13.09.12  | End time:           | 14:39:30 |

#### Titration data

|                    |                     |                  |                       |
|--------------------|---------------------|------------------|-----------------------|
| Sample ID:         | TRIS                | Weight:          | 0.1038 g              |
| Start pH:          | pH 9.590            | End pH:          | pH 3.864              |
| Start temperature: | 25.0 °C (m)         | End temperature: | 25.0 °C (m)           |
| Zero point:        | pH 6.83 / -10.0 mV  | Slope:           | 100.6 % / -59.5 mV/pH |
| EQ:                | 8.560 ml / pH 4.849 | Titre:           | 0.1001 mol/l          |

#### Calculation formula

Titre:  $(W \cdot F2) / ((EQ1 - B) \cdot M \cdot F1) \rightarrow M103$   
 Mol (M): 121.14000

|                  |           |                |           |
|------------------|-----------|----------------|-----------|
| Weight (W):      | man       | Factor 2 (F2): | 1000.0000 |
| Blank value (B): | 0.0000 ml | Factor 1 (F1): | 1.0000    |
| Statistics:      | Off       |                |           |

# Application

Page 2: Method parameters Titer determination:

## Method data overall view

|                          |                     |                       |                   |
|--------------------------|---------------------|-----------------------|-------------------|
| Method name:             | Titre HCl           | Created at:           | 09/13/12 14:23:02 |
| Method type:             | Automatic titration | Last modification:    | 09/13/12 14:27:56 |
| Measured value:          | pH                  | Damping settings:     | None              |
| Titration mode:          | Dynamic             | Documentation:        | GLP               |
| Dynamic:                 | Steep               |                       |                   |
| Measuring speed / drift: | Normal:             | minimum holding time: | 02 s              |
|                          |                     | maximum holding time: | 15 s              |
|                          |                     | Measuring time:       | 02 s              |
|                          |                     | Drift:                | 20 mV/min         |
| Initial waiting time:    | 0 s                 |                       |                   |
| Titration direction:     | Decrease            |                       |                   |
| Pretitration:            | Off                 |                       |                   |
| End value:               | 2.500 pH            |                       |                   |
| EQ:                      | On (1)              |                       |                   |
| Slope value:             | Steep               | Value:                | 700               |

## Dosing parameter

|                        |          |                |      |
|------------------------|----------|----------------|------|
| Dosing speed:          | 100 %    | Filling speed: | 30 s |
| Maximum dosing volume: | 50.00 ml |                |      |

## Unit values

|                          |                   |
|--------------------------|-------------------|
| Unit size:               | 20ml              |
| Unit ID:                 | 10039005          |
| Reagent:                 | HCl 0.1 mol/L     |
| Batch ID:                | no Charge         |
| Concentration [mol/l]:   | 0.10070           |
| Determined at:           | 12/05/11 19:18:45 |
| Expire date:             | 08/18/12          |
| Opened/compounded:       | 09/10/11          |
| Test according ISO 8655: | 05/10/11          |
| Last modification:       | 09/13/12 14:35:18 |

# Application

---

## Di n-butylamine solution

129 g di-n-butyl amine are weight in a 1 L volumetric flask and filled up to the mark with dried toluene.

## Titration of the sample

There are different procedures. Sometimes the di n-butylamine solution is pipetted first and then the sample is weight in. In some procedures the sample is weight in first and the reagent mixture is added afterwards.

The ISO 14896 (part A) describes the pipetting of the di n-butylamine solution into an Erlenmeyer flask (250-300 ml) and then the addition of the sample. We are following this procedure.

Pipette 50.00 ml using a volumetric pipette into an Erlenmeyer flask with stopper. The sample is weight in using a 2 or 5 ml syringe by different weight. Sample weight should be between 2- 5 g.

Close the flask with a stopper and gently stir the mixture on the magnetic stirrer until the sample is completely dissolved. Wait then additional 15 min for the reaction + 5 – 10 minutes until the mixture cooled down to room temperature. Add 100 – 150 ml acetone (or any other suitable solvent) carefully.

Place the Erlenmeyer flask on the magnetic stirrer (if not already) and start the titration method. After the titration rinse the electrode and burette tip with ethanol. For each set of samples perform a blank titration.

## Result calculation

% NCO

$$\% \text{ NCO} = (EQ1-B) * M * T * F1/(W*F2)$$

EQ1: ml consumption at the equivalence point

B: ml consumption for the blank titration

M: 4.202 equivalent weight including conversation factor to % NCO

T: concentration of the HCl titrant (e.g. 1.002 mol/l)

W: sample weight in g


F1,F2: 1

## Application

Blank titration page 1: Curve and result

### GLP documentation

#### Titration graph



#### Method data

|              |                     |                     |          |
|--------------|---------------------|---------------------|----------|
| Method name: | NCO blank titration | Titration duration: | 2 m 13 s |
| End date:    | 27.06.12            | End time:           | 10:53:46 |

#### Titration data

|            |                      |         |           |
|------------|----------------------|---------|-----------|
| Sample ID: | Test 2               | Weight: | 1.0000 g  |
| Start mV:  | -173.3 mV            | End mV: | 272.0 mV  |
| EQ:        | 16.330 ml / 151.2 mV | Blank:  | 16.330 ml |

#### Calculation formula

|          |            |
|----------|------------|
| Blank:   | EQ1 -> M01 |
| Mol (M): | 1.00000    |

Statistics: Off

# Application

Blank titration page 2: method

The pre-titration volume can also be changed to 25-40 ml. The end volume can also be changed to 55 or 60 ml.

## Method data overall view

|                          |                     |                       |                   |
|--------------------------|---------------------|-----------------------|-------------------|
| Method name:             | NCO blank titration | Created at:           | 11/17/11 11:29:35 |
| Method type:             | Automatic titration | Last modification:    | 06/27/12 10:48:12 |
| Measured value:          | mV                  |                       |                   |
| Titration mode:          | Dynamic             | Documentation:        | GLP               |
| Dynamic:                 | Steep               |                       |                   |
| Measuring speed / drift: | Normal:             | minimum holding time: | 02 s              |
|                          |                     | maximum holding time: | 15 s              |
|                          |                     | Measuring time:       | 02 s              |
|                          |                     | Drift:                | 20 mV/min         |
| Initial waiting time:    | 0 s                 |                       |                   |
| Titration direction:     | Increase            |                       |                   |
| Pretitration:            | 10.000 ml           | Delay time:           | 10 s              |
| End value:               | Off                 |                       |                   |
| EQ:                      | On                  |                       |                   |
| Slope value:             | Steep               | Value:                | 700               |

## Dosing parameter

|                        |          |                |      |
|------------------------|----------|----------------|------|
| Dosing speed:          | 100 %    | Filling speed: | 30 s |
| Maximum dosing volume: | 50.00 ml |                |      |

## Unit values

|                          |                   |
|--------------------------|-------------------|
| Unit size:               | 20ml              |
| Unit ID:                 | 10039060          |
| Reagent:                 | HCl 1 mol/l       |
| Batch ID:                | no entry          |
| Concentration [mol/l]:   | 1.00000           |
| Determined at:           | 06/27/12 16:03:58 |
| Expire date:             | --                |
| Opened/compounded:       | --                |
| Test according ISO 8655: | --                |
| Last modification:       | 06/27/12 9:04:36  |


## Application

Sample titration page 1: Curve and result

The result calculation was wrong here for this titration because of using a wrong factor.

### GLP documentation

#### Titration graph



#### Method data

Method name: NCO sample  
End date: 27.06.12

Titration duration: 2 m 11 s  
End time: 11:03:06

#### Titration data

Sample ID: Dummy one  
Start mV: -133.7 mV

Weight: 5.5463 g  
End mV: 271.1 mV

EQ: 10.448 ml / 155.1 mV

free NCO: 942.928 %



# Application

Sample titration page 2: method

The pre-titration volume can also be changed to 5, 10 or more ml.

## Method data overall view

|                          |                     |                       |                   |
|--------------------------|---------------------|-----------------------|-------------------|
| Method name:             | NCO sample          | Created at:           | 12/06/11 12:43:30 |
| Method type:             | Automatic titration | Last modification:    | 06/27/12 10:24:02 |
| Measured value:          | mV                  |                       |                   |
| Titration mode:          | Dynamic             | Documentation:        | GLP               |
| Dynamic:                 | Steep               |                       |                   |
| Measuring speed / drift: | Normal:             | minimum holding time: | 02 s              |
|                          |                     | maximum holding time: | 15 s              |
|                          |                     | Measuring time:       | 02 s              |
|                          |                     | Drift:                | 20 mV/min         |
| Initial waiting time:    | 0 s                 |                       |                   |
| Titration direction:     | Increase            |                       |                   |
| Pretitration:            | 2.000 ml            | Delay time:           | 10 s              |
| End value:               | Off                 |                       |                   |
| EQ:                      | On                  |                       |                   |
| Slope value:             | Steep               | Value:                | 700               |

## Dosing parameter

|                        |          |                |      |
|------------------------|----------|----------------|------|
| Dosing speed:          | 100 %    | Filling speed: | 30 s |
| Maximum dosing volume: | 50.00 ml |                |      |

## Unit values

|                          |                   |
|--------------------------|-------------------|
| Unit size:               | 20ml              |
| Unit ID:                 | 10039060          |
| Reagent:                 | HCl 1 mol/l       |
| Batch ID:                | no entry          |
| Concentration [mol/l]:   | 1.00000           |
| Determined at:           | 06/27/12 16:03:58 |
| Expire date:             | --                |
| Opened/compounded:       | --                |
| Test according ISO 8655: | --                |
| Last modification:       | 06/27/12 9:04:36  |

## Application

---

### Notes

If you have any questions on the application, you can feel free to contact us..

SI Analytics GmbH  
Hattenbergstr. 10  
55122 Mainz  
Germany

Phone: +49 (0) 6131 / 66 – 5062  
+49 (0) 6131 / 66 – 5118  
Fax: +49 (0) 6131 / 66 – 5001  
E-Mail: [titration@si-analytics.com](mailto:titration@si-analytics.com)  
Homepage: [www.si-analytics.com](http://www.si-analytics.com)