

## Application

---

**Determination of acidity/acid  
degree in milk and milk  
products**

## Application

---

### Use

This method is used for the quantitative determination of acidity (acid degree) in milk and milk products. The acidity is calculated either as SH-, Dornic- or Therner degree. The different degree uses different concentrations of NaOH titrant.

### Appliances

- Titrator: TL 6000/7000 (TL 6000/7000 M2/20) consists of
- Basic device
- Magnetic stirrer TM 235
- 20 mL Exchange unit WA 20, with brown glass bottle for titrant complete
- And pH combination electrode A 162 DIN ID

### Electrodes

- Electrode: A 162 DIN ID
- Calibration: DIN buffer pH= 4.00 and pH= 7.00

# Application

## Reagents

- Titrant: sodium hydroxide solution 0.1 mol/l (Therner), 0.25 (SH) or 0.11 (Dornic) mol/l
- Soda lime for carbon dioxide uptake of the reagent.
- Titer: potassium hydrogen phthalate (reference material)

## Description

### Calibration

The pH combination electrode is calibrated in technical buffer pH=4.00 and pH= 7.00 or in DIN buffer pH= 4.01 and pH= 6.87.

Example of the calibration documentation:

### Calibration

#### Buffers used

pH buffer 1:	TEC_4.000
pH buffer 2:	TEC_7.000

#### Measured values

pH buffer 1:	TEC_4.000	165.6 mV / 23.4 °C
pH buffer 2:	TEC_7.000	-11.2 mV / 23.0 °C

#### Calibration data

Slope:	99.4 % / -58.8 mV/pH
Zero point:	pH 6.81 / -11.2 mV
Temperature:	23.4 °C (a)
Date and time:	07.03.13 / 15:04

### Determination of the exact concentration of the standard solution

By carbon dioxide absorption from the air occurs in the sodium hydroxide solution of sodium bicarbonate, which changes the pH of the titrant. To prevent this, a drying tube filled with soda lime is placed on the reagent bottle. The exact concentration of the sodium hydroxide solution is determined using the standard potassium hydrogen phthalate. The potassium hydrogen phthalate is dried in the oven before the titer determination for 2 hours at 120°C and cooled in a desiccator.

#### Implementation


In a 50 mL beaker, 0.2 to 0.3g potassium hydrogen phthalate were weighed accurately and dissolved in 30 mL of dist. water with stirring. It is titrated with 0.1 mol/l sodium hydroxide solution. For 0.25 mol/l NaOH you need 0.5 – 0.75 g potassium hydrogen phthalate.

## Application

### Standardisation titration (result)

#### GLP documentation

##### Titration graph



##### Method data

Method name:	Titre NaOH	Titration duration:	2 m 15 s
End date:	08.01.13	End time:	15:46:03

##### Titration data

Start pH:	pH 4.065	Weight:	0.20490 g
Start temperature:	25.0 °C (m)	End pH:	pH 9.667
Zero point:	pH 6.85 / -8.9 mV	End temperature:	25.0 °C (m)
EQ:	10.032 ml / pH 8.498	Slope:	98.7 % / -58.4 mV/pH
Mean value:	---	Titre:	0.1000 mol/l
		RSD:	---

##### Calculation formula

Titre:	$(W \cdot F2) / ((EQ1 - B) \cdot M \cdot F1) \rightarrow WA$	Mol (M):	204.22000
--------	--	----------	-----------

Weight (W):	0.2049 g (m)	Factor 2 (F2):	1000.0000
Blank value (B):	0.0000 ml	Factor 1 (F1):	1.0000
Statistics:	3		

## Application

### Standardisation titration (method)

#### Method data

Method name:	Titer NaOH	Created at:	02/15/12 15:32:03
Method type:	Automatic titration	Last modification:	02/16/12 10:14:55
Measured value:	pH	Documentation:	GLP
Titration mode:	Dynamic		
Dynamic:	average		
Measuring speed / drift:	Normal:	minimum holding time:	02 s
		maximum holding time:	15 s
		measuring time:	02 s
		drift:	20 mV/min
Initial waiting time:	0 s		
Titration direction:	Increase		
Pretitration:	Off		
End value:	10.500 pH		
EQ:	On		
slope value:	Steep	Value:	700

#### Dosing parameter

Dosing speed:	100 %	Filling speed:	30 s
Maximum dosing volume:	30.00 ml		

#### Calculation formula

Titer NaOH 0,1mol/l:	$(W * F2) / ((EQ1 - B) * M * F1)$		
Mol (M):	20.42230		
Unit:		Decimal places:	4

Weight (W):	man	Factor 2 (F2):	1000.0000
Blank value (B):	0.0000 ml	Factor 1 (F1):	1.0000

#### Device information

Device: TitroLine 6000

Serial number:

Software version: 07\_12

mth\_Titer\_NaOH\_29\_02\_12-10\_51\_09.pdf

1/1

## Application

### Titration of the sample

Pipette accurately in 100 ml beaker either:


- 25 ml milk for SH degree
- 11 ml milk for Dornic degree
- 9 ml for Therner degree

Add app. 25 ml – 40 ml CO<sub>2</sub> free water into and mix the sample on the magnetic stirrer for a few seconds. Titrate with the suitable NaOH titrant to a fixed pH endpoint (8.2-8.7, depending on national norms).

### Example

#### GLP documentation

##### Titration graph



##### Method data

Method name:	acidity in milk	Titration duration:	2 m 5 s
End date:	29.04.13	End time:	13:47:45

##### Titration data

Sample ID:	Milk 4	Weight:	24.67000 g
Start pH:	pH 6.768	End pH:	pH 8.215
Start temperature:	19.1 °C (a)	End temperature:	19.3 °C (a)
Zero point:	pH 6.82 / -10.3 mV	Slope:	98.8 % / -58.4 mV/pH
EP1:	1.806 ml / pH 8.200	Acidity:	7.321 °SH

##### Calculation formula

Acidity:	$(EP1-B)*T*M*(F1)/(W*F2)$	Mol (M):	1.00000
Blank value (B):	0.0000 ml	Titre (T):	1.00000000 (a)
Factor 1 (F1):	100.0000	Weight (W):	24.67000 g (m)
Factor 2 (F2):	1.0000	Statistics:	Off

# Application

## Method parameters: (optimized for SH degree)

### Method data overall view

Method name:	acidity in milk	Created at:	04/29/13 13:24:34
Method type:	Automatic titration	Last modification:	04/29/13 13:44:18
Measured value:	pH	Damping settings:	None
Titration mode:	End pt.	Documentation:	GLP
Linear steps:	0.050 ml		

Measuring speed / drift:	Normal:	minimum holding time:	02 s
		maximum holding time:	15 s
		Measuring time:	02 s
		Drift:	20 mV/min
Initial waiting time:	0 s		
Titration direction:	Increase		
Pretitration:	Off		
Endpoint 1:	pH 8.200	delta endpoint 1:	pH 1.000
		Endpoint delay 1:	5 s
Endpoint 2:	Off		

### Dosing parameter

Dosing speed:	20.00 %	Filling speed:	30 s
Maximum dosing volume:	50.00 ml		

### Unit values

Unit size:	20ml
Unit ID:	00360005
Reagent:	NaOH
Batch ID:	no Charge
Concentration [mol/l]:	1.00000
Determined at:	04/29/13 20:26:19
Expire date:	01/01/13
Opened/compounded:	02/01/12
Test according ISO 8655:	01/01/00
Last modification:	04/29/13 13:26:21

### Device information

Device:	TitroLine 7750
Serial number:	10018602
Software version:	1316

Hard Disk\result\acidity\_in\_milk\_29\_04\_13-13\_45\_40.pdf

2/2

# Application

---

## Notes

The acidity can also calculate as % lactic acid instead of acid degree. In this case the calculation is:

% lactic acid:  $(EP1-B) * M * T * F1 / (W * F2)$

EP1: ml consumption to pH endpoint

B: blank value in ml. Here = 0

M: molecular weight of lactic acid: 90.08

T: exact concentration of the NaOH in mol/l (e.g. 0.2510 mol/l)

F1: 0.1 (conversion factor to %)

W: sample amount in g

F2: 1

If you have any questions on the application, you can feel free to contact us.

SI Analytics GmbH  
Hattenbergstr. 10  
55122 Mainz  
Germany

Phone: +49 (0) 6131 / 66 – 5062  
+49 (0) 6131 / 66 – 5118  
Fax: +49 (0) 6131 / 66 – 5001  
E-Mail: [titration@si-analytics.com](mailto:titration@si-analytics.com)  
Homepage: [www.si-analytics.com](http://www.si-analytics.com)