Freezer Flask

Purpose

To illustrate an endothermic reaction involving reactants in the solid state.

Materials

- Barium hydroxide 8-hydrate, Ba(OH)₂·8H₂O
- Ammonium thiocyanate, NH₄SCN
- 500-mL Erlenmeyer flask with rubber stopper
- Wood board
- Electronic balance
- Water

Safety

- Read the SDS sheets for all chemicals before using them.
- Ammonia gas is produced. Avoid inhalation. (You may wish to do this demonstration in a fume hood)
- Wear safety glasses, gloves, and lab coat.
- Firmly hold the rubber stopper to prevent it from popping off.

Procedure

- 1. Pour a small amount of water (enough to ensure good contact with the base of the Erlenmeyer flask) on a square wooden board.
- 2. Place 50 g of barium hydroxide octahydrate and 25 g ammonium thiocyanate into the Erlenmeyer flask. Place a rubber stopper in the flask. Swirl to mix.
- 3. Place the flask in the puddle of water on the board.
- 4. Lift the flask after a few minutes.

Results

• The reaction vessel gets cold enough to freeze the water between the board and the flask.

Follow-up Teaching Notes

- Energy absorbed by the reaction freezes the water so that the board will be attached to the flask.
- Reaction is

 $Ba(OH)_2 \cdot 8H_2O_{(s)} + 2 NH_4SCN_{(s)} \rightarrow Ba(SCN)_{2(aq)} + 2 NH_{3(q)} + 10 H_2O_{(l)}$

• Ammonium nitrate can be a substitute for ammonium thiocyanate.

Connections

• Enthalpy and heat of reaction, reaction between two solids.

Extensions

- Demonstration of an endothermic reaction: dissolve ammonium chloride in water.
- Demonstration of an exothermic reaction: dissolve calcium chloride in water.

Disposal/Clean-up

- Open the flask in a fume hood and allow the ammonia to evaporate overnight.
- Remaining contents can be placed in the science department's heavy metal waste container for proper disposal.

Ward's Science Tel: (866) 260-0501