# **Under Pressure**

#### **Purpose**

To illustrate expansion of gases under reduced pressure.

#### **Materials**

- Shaving cream
- Marshmallow
- Soda water
- Doll's head
- Balloon
- Tape
- Vacuum pump with bell jar

## Safety

- Wear safety glasses and gloves.
- Ensure there are no cracks in the bell jar.

#### **Procedure**

- 1. Fill an empty doll's head with shaving cream from the neck opening. Ensure openings are present at the eyes and ears.
- 2. Tape the neck opening closed.
- 3. Place the doll's head in a bell jar connected to vacuum pump.
- 4. Evacuate the chamber of air.
- 5. Other items that can be placed in the bell jar include soda water, small balloon, or a marshmallow.

#### Results

• As the surrounding pressure decreases, the shaving cream/marshmallow expands.

### **Follow-up Teaching Notes**

- Shaving cream is a foam (type of colloid) that consists of a gas evenly distributed in a liquid.
- As the surrounding pressure is decreased, the gas expands in size.

#### **Connections**

Gas laws, solubility of gases, decompression sickness, colloids

#### Extension

Can be related to decompression sickness.

#### Disposal/Clean-up

- The marshmallow and shaving cream can be disposed of in the garbage.
- The soda water can be washed down the drain.

Ward's Science Tel: (866) 260-0501

114